

ZrGEN
Ti-Base



Library info		Ürün Bilgisi				
ZrGEN Option In exocad library		Çap (Ø)	Cuff Yükseklik (mm)	Post Yükseklik (mm)	Tip	Ref.C
System	Kit					
	ZrGEN 4015	Ø4.0	0.6	4.5	Hex	AAOIPR4015.MTN
	ZrGEN 4016			6.0		AAOIPR4016.MTN
	ZrGEN 4018			8.0		AAOIPR4018.MTN
	ZrGEN 4025		1.5	4.5		AAOIPR4025.MTN
	ZrGEN 4026			6.0		AAOIPR4026.MTN
	ZrGEN 4028			8.0		AAOIPR4028.MTN
	ZrGEN 4035		3.0	4.5		AAOIPR4035.MTN
	ZrGEN 4036			6.0		AAOIPR4036.MTN
	ZrGEN 4038			8.0		AAOIPR4038.MTN
	ZrGEN 4045		4.0	4.5		AAOIPR4045.MTN
	ZrGEN 4046			6.0		AAOIPR4046.MTN
	ZrGEN 4048			8.0		AAOIPR4048.MTN
	ZrGEN 4015	Ø4.0	0.6	4.5	Non-Hex	AAOIPR4015N.MTN
	ZrGEN 4016			6.0		AAOIPR4016N.MTN
	ZrGEN 4018			8.0		AAOIPR4018N.MTN
	ZrGEN 4025		1.5	4.5		AAOIPR4025N.MTN
	ZrGEN 4026			6.0		AAOIPR4026N.MTN
	ZrGEN 4028			8.0		AAOIPR4028N.MTN
	ZrGEN 4035		3.0	4.5		AAOIPR4035N.MTN
	ZrGEN 4036			6.0		AAOIPR4036N.MTN
	ZrGEN 4038			8.0		AAOIPR4038N.MTN
	ZrGEN 4045		4.0	4.5		AAOIPR4045N.MTN
	ZrGEN 4046			6.0		AAOIPR4046N.MTN
	ZrGEN 4048			8.0		AAOIPR4048N.MTN
	ZrGEN 4515	Ø4.5	0.6	4.5	Hex	AAOIPR4515.MTN
	ZrGEN 4516			6.0		AAOIPR4516.MTN
	ZrGEN 4518			8.0		AAOIPR4518.MTN
	ZrGEN 4525		1.5	4.5		AAOIPR4525.MTN
	ZrGEN 4526			6.0		AAOIPR4526.MTN
	ZrGEN 4528			8.0		AAOIPR4528.MTN
	ZrGEN 4535		3.0	4.5		AAOIPR4535.MTN
	ZrGEN 4536			6.0		AAOIPR4536.MTN
	ZrGEN 4538			8.0		AAOIPR4538.MTN
	ZrGEN 4545		4.0	4.5		AAOIPR4545.MTN
	ZrGEN 4546			6.0		AAOIPR4546.MTN
	ZrGEN 4548			8.0		AAOIPR4548.MTN
	ZrGEN 4515	Ø4.5	0.6	4.5	Non-Hex	AAOIPR4515N.MTN
	ZrGEN 4516			6.0		AAOIPR4516N.MTN
	ZrGEN 4518			8.0		AAOIPR4518N.MTN
	ZrGEN 4525		1.5	4.5		AAOIPR4525N.MTN
	ZrGEN 4526			6.0		AAOIPR4526N.MTN
	ZrGEN 4528			8.0		AAOIPR4528N.MTN
	ZrGEN 4535		3.0	4.5		AAOIPR4535N.MTN
	ZrGEN 4536			6.0		AAOIPR4536N.MTN
	ZrGEN 4538			8.0		AAOIPR4538N.MTN
	ZrGEN 4545		4.0	4.5		AAOIPR4545N.MTN
	ZrGEN 4546			6.0		AAOIPR4546N.MTN
	ZrGEN 4548			8.0		AAOIPR4548N.MTN



* ZrGEN Library provides Cement gap (ZrGEN and Abutment space) 0, 10, 20, 30, 40, 50 µ.